



EX-RNA™-SPINTUBE

RNA Extraction Kit For RT-PCR Testing (Spintube Method) (For *in vitro* diagnostic use only)

INTENDED USE

This kit is recommended to be used for the isolation of viral RNA from various samples like fresh and frozen plasma, serum, nasopharyngeal swab, oropharyngeal swab, sputum, broncho-alveolar lavage in viral transport medium and other body fluids. The system is intended for use by professional users trained in molecular biological techniques.

SUMMARY

EX-RNA™-SPINTUBE is a simple and efficient system for purification of total RNA from samples like fresh and frozen plasma, serum, nasopharyngeal swab, oropharyngeal swab, sputum, broncho-alveolar lavage in viral transport medium and other body fluids. Ribonucleic acid (RNA) purified from biological materials is utilized extensively for molecular biology research and is becoming an important tool in human clinical testing.

PRINCIPLE

Samples are lysed and homogenized in lysis buffer, which contains chaotropic salts. Chaotropic salts have two important roles in nucleic acid extraction. Firstly, they destabilize hydrogen bonds, van der Waals forces and hydrophobic interactions, leading to destabilization of proteins, including nucleases. Secondly, they disrupt the association of nucleic acids with water, thereby providing optimal conditions for binding. Buffering conditions enhances the binding of nucleic acids in the column. After centrifuging the lysate through the silica membrane, the desired nucleic acids are bound to the column and unbound components are in the flow-through. However, the membrane will contain impurities such as protein, polysaccharides and salt residues. The wash steps remove such impurities. Subsequent addition of elution buffer to the column will hydrate the nucleic acids which is collected in the elute collection tube.

PRESENTATION

REF	1108100100
▽	100 Tests

Kit Components (for 100 Tests):

Reagents:

R1: Lysis Buffer 60 mL
R2: Wash Buffer 30 mL
R3: Elution Buffer 10 mL

Accessories:

T1: Sample Lysis Tubes 100 Nos.
T2: Spin Column with Wash Tubes 100 Nos.
T3: Elute Collection Tubes 100 Nos.

Packinsert 1 No.

STORAGE AND STABILITY OF THE KIT

The kit and reagents are to be stored at 20-25°C.

REAGENT/MATERIALS REQUIRED BUT NOT PROVIDED

- Ethanol (96-100%)



- RNase inactivating agents such as **DECON-R™** (10 x 200 mL - Cat. No. 1108092000)
- Micro-centrifuge (operational requirement 2500 x g - 6000 x g)
- Vortex
- Microtube Stand
- Sterile Nuclease free microtips
- Micropipettes

WARNINGS

EX-RNA™-SPINTUBE contains guanidine thiocyanate which is corrosive to metals, causes skin corrosion and serious eye damage.

- Recommended Personal Protective Equipment Includes Dust mask type N95, Eye shields and thick durable Nitrile or Plastic Gloves.
- If on skin: Gently wash with plenty of water.
- If skin irritation or rash occurs: Get medical advice/attention.

Kindly note that this transport media should not be used in a testing platform that uses bleach or in laboratories that use bleach as a part of their routine decontamination and disposal process. When the bleach interacts with guanidine thiocyanate it produces the hazardous cyanide gas.

PRECAUTIONS

- A thorough understanding of the pack insert is mandatory before performing the test for the first time. Adherence to protocol specified herein is necessary to ensure optimal performance of the product. Any deviation from the assay procedure may affect the results.
- Do not mix reagents of different lots.
- Do not use reagents beyond their expiry date.
- Use nuclease free plasticware and water.
- Avoid any contamination among samples; for this purpose, disposable tips should be used for each sample and reagent.
- Do not use reagents from other manufacturers along with the kit reagents for a given test run.
- Do not interchange reagent vials and their screw caps to avoid cross contamination. Use a clean, fresh, disposable pipette tip for each reagent or specimen manipulation.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- Always wear gloves. After wearing gloves, do not touch surfaces and equipment to avoid reintroduction of RNases to material.
- Treat surfaces of benches and glassware with RNase inactivating agents such as **DECON-R™** (10 x 200 ml - Cat No. 1108092000).
- Use sterile, disposable RNase-free plastic ware.
- Use only DEPC-treated water or nuclease free water such as **RNase FREE WATER™** (10 x 125 ml - Cat No. 1108092000, 10 x 500ml - Cat. No. 1108115000).
- Best results were obtained by using fresh samples or samples that have been quickly frozen in liquid nitrogen or stored at -70°C.
- Practice Good Laboratory practice (GLP) while handling specimen and other reagents.
- Follow precautions as per the handling of specimen

- capable of transmitting infectious agents.
- Lysis buffer contains chaotropic salt. It is not compatible with disinfecting agents that contain bleach. Dispose it as per local medical waste guidelines.

PREPARATION OF REAGENTS:

Reagents provided as a concentrate

Reagents	Contents/Labels	Volume of Ethanol (96-100 %) to be added
R2	Wash Buffer	30 mL Ethanol

Reagents Not Provided, to be prepared by the User

Reagents	Contents/Labels	Volume of Ethanol (96-100 %) to be added
B1	Binding Reagent	60 mL Ethanol
W2	Ethanol Wash Solution	45 mL Ethanol + 15 mL Nuclease free water

Note: R2 is supplied as concentrate. B1 and W2 to be prepared separately by user. Before using for the first time, add recommended volume of ethanol (96–100%) as indicated in the table above and also label them on the respective bottles. Please tick mark the check box provided on the label after addition of ethanol.

PROCEDURE:

- Take 500 µL of **R1, Lysis Buffer** in the sample lysis tube and add 300 µL of sample to it. Then add 500 µL of **B1, Binding Reagent**.
Optional step - Poly(A) RNA: Pipette 10 µL Poly(A) RNA and add it to the reaction mix prepared in step 1.
- Mix well the reaction mix by pipetting up and down 10 to 15 times and incubate for 10 minutes at room temperature.
- Transfer 650 µL of this suspension into the spin column fitted in wash tube.
- Centrifuge at 6000 x g for 1 minute at room temperature. Discard the flow-through from wash tube and blot dry on absorbent paper.
- Add the remaining 650 µL suspension from the sample lysis tube into the same spin column fitted in wash tube.
- Centrifuge at 6000 x g for 1 minute at room temperature. Discard the flow-through from wash tube and blot dry on absorbent paper.

Symbol keys

25°C 20°C Store at 20-25°C	Manufacturer	In vitro Diagnostic Medical Device	R1 Lysis Buffer	T1 Sample Lysis Tubes	Acute Toxicity
Use by (Last day of stated month)	Consult Instructions for use	Batch Number	R2 Wash Buffer(4x)	T2 Spin Column with Wash Tubes	Health hazard
Date of Manufacture	Catalogue Number	Authorised Representative in the European Community	R3 Elution Buffer	T3 Elute Collection Tubes	This way up Corrosive

0720/VER-01



Manufactured by:

Coral Clinical Systems

A Division of Tulip Diagnostics (P) Ltd.

PLOT NO. M-46, BLDG. NO. D, PHASE III B, VERNA IND. EST., VERNA, GOA-403 722, INDIA.
REGD. OFFICE : GITANJALI, TULIP BLOCK, DR. ANTONIO DO REGO BAGH,
ALTO SANTACRUZ, BAMBOLIM COMPLEX P.O., GOA-403 202, INDIA.

CMC Medical Devices & Drugs S.L.,
C/ Horacio Lengo No. 18, CP 29006, Malaga, Spain.

EC REP