

Hayem's (RBC Diluting) Fluid

Intended Use

Hayem's (RBC Diluting) Fluid is used for RBC count.

Summary & Principle

RBC diluting fluid is isotonic with blood; hence hemolysis does not take place. Normal Saline also can be used. But it causes slight creation of red blood cells and allows rouleaux formation. The blood specimen is diluted 1:200 with the RBC diluting fluid and cells are counted under high power (40x objective) by using a counting chamber. The number of cells in undiluted blood are calculated and reported as the number of red cells per cu mm (MI) of whole blood.

Reagent/ Contents

Mercuric chloride	0.25 gm
Sodium sulphate	2.50 gm
Sodium chloride	0.50 gm
Distilled water	100.0mL

**Formula adjusted, standardized to suit performance parameters

Directions

1. Draw EDTA anticoagulated blood to exactly the 0.5 mark of the RBC pipette.
2. Wipe the tip of the pipette, clean with a piece of dry gauze without touching the opening of the capillary and immerse in the freshly filtered diluting fluid.
3. Do not insert the pipette in the bottle of counting solution.
4. By gentle mouth suction, draw the diluting fluid steadily into the pipette to exactly the 101 mark past the bulb, rotating the pipette on its long axis to ensure thorough mixing of blood and diluent.
5. Immediately mix the contents of the pipette thoroughly by placing the thumb over one end and shake for 1 minute.
6. Diluted blood must be examined within 2 hours.

Appearance

Colourless, solution.

Storage and Stability

Store at 15°C-30°C away from bright light. Stability of the Hayem's (RBC Diluting) Fluid is as per expiry date mentioned on label.

Warranty

This product is designed to perform as described on the label and package insert. The manufacturer disclaims any implied warranty of use and sale for any other purpose.

Reference

1. Data on file: Microxpress[®], A Division of Tulip Diagnostics (P) Ltd.

Product Presentation

Cat. No.	Product Description	Pack Size
207080090500	Hayem's (RBC Diluting) Fluid	500 mL

Disclaimer

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.
