



# Clinical Microbiology Products

## Selective Culture Media Recommended for Bacteria

### **Bacillus**

- Bacillus Cereus Agar Base
- MYP Agar

### **Clostridium**

- Columbia Blood Agar Base
- Cooked Meat Medium
- Reinforced Clostridial Agar
- Reinforced Medium for Clostridia

### **Corynebacterium**

- Columbia Blood Agar Base
- Hartley's Digest Broth
- Loeffler Medium Base

### **Legionella**

- Buffered Charcoal Yeast Extract Agar
- Legionella Agar Base

### **Klebsiella**

- C.L.E.D. Agar
- Michrom™ UTI Agar
- Decarboxylase Agar Base
- EMB Agar
- Endo Agar
- MacConkey Agar with Crystal Violet, NaCl and 0.15% Bile Salts

### **Enterobacteriaceae (Enrichment)**

- Brilliant Green Bile Lactose Broth 2%
- Buffered Peptone Water
- E.C. Broth
- Enterobacteria Enrichment Broth, Mossel
- Lactose Broth

### **Enterobacteriaceae (Isolation)**

- C.L.E.D. Agar
- EMB Agar
- Michrom™ Coliform Agar
- Michrom™ E.coli Agar
- Michrom™ UTI Agar
- MacConkey Agar with Crystal Violet, NaCl and 0.15% Bile Salts
- MacConkey Agar without Crystal Violet and with 0.15% Bile Salts
- Simmons Citrate Agar
- Violet Red Bile Agar

### **Mycobacterium**

- Lowenstein Jensen Medium Base
- Middlebrook 7H9 Agar Base
- Middlebrook 7H10 Agar Base
- Middlebrook 7H11 Agar Base

### **Pseudomonas**

- Cetrimide Agar Base
- King's Medium A Base
- King's Medium B Base
- Pseudomonas Agar Base
- Pseudomonas Agar Medium for Detection of Fluorescein USP
- Pseudomonas Agar Medium for Detection of Pyocyanin USP

### **Proteus**

- C.L.E.D. Agar
- Michrom UTI Agar
- Phenylalanine Agar
- Urea Agar Base, Christensen
- Urea Broth Base

### **Staphylococcus**

- Baird Parker Agar Base
- Dnase Test Agar Base
- Mannitol Salt Agar
- Salt Meat Broth
- Tellurite Glycine Agar Base
- Vogel Johnson Agar Base w/o Tellurite

### **Streptococcus**

- Bile Esculin Agar
- Bile Esculin Azide Agar
- Blood Agar Base (Infusion Agar)
- Brain Heart Infusion Agar
- Brain Heart Infusion Broth
- Columbia Blood Agar Base
- Todd Hewitt Broth

### **Salmonella (Enrichment)**

- Fluid Selenite Cystine Medium (Twin Pack)
- Mueller Kauffmann Tetrathionate Broth Base (Novobiocin Broth)
- Rappaport Vassiliadis Salmonella Enrichment Broth
- Selenite F Broth (Twin Pack)
- Tetrathionate Brilliant Green Bile Broth
- Tetrathionate Broth Base, Hajna

### **Salmonella (Isolation)**

- Bismuth Sulphite Agar
- Brilliant Green Agar, Modified
- C.L.E.D. Agar
- Deoxycholate Citrate Agar
- Endo Agar
- Hektoen Enteric Agar
- Lysine Iron Agar
- SS (Salmonella Shigella Agar) Agar
- Triple Sugar Iron Agar
- Wilson Blair Agar Base
- XLD Agar

### **Shigella**

- GN Broth
- SS (Salmonella Shigella Agar) Agar

### **Vibrio**

- Gelatin Agar
- TCBS Agar

**Citrate Utilization Test**

Simmons Citrate Agar  
Yeast Nitrogen Agar  
Yeast Nitrogen Base Agar (Twin Pack)

**Esculin Hydrolysis Test**

Bile Esculin Agar  
Bile Esculin Azide Agar

**Hemolysis Test**

Blood Agar Base (Infusion Agar)  
Brain Heart Infusion Agar  
Columbia Blood Agar Base  
Nutrient Agar  
Soyabean Casein Digest Agar



**H<sub>2</sub>S Production Test**

Kligler Iron Agar  
Lysine Iron Agar  
SIM Medium  
Triple Sugar Iron Agar



**Indole Production Test**

Lauryl Tryptose Mannitol Broth w/ Tryptophan  
Peptone Water  
Tryptone Water (Tryptone Broth)

**Lecithinase Activity**

Bacillus Cereus Agar Base  
Baird Parker Agar Base  
Mannitol Salt Agar

**Motility Test**

Mannitol Motility Test Medium  
Motility Test Medium

**Urease Test**

Urea Agar Base, Christensen  
Urea Broth Base

**Biochemical Identification Kits**

- Bifidobacterium Identification Kit
- Candida Identification Kit
- E. coli Identification Kit
- Enterobacteriaceae Identification Test Kit
- IMViC Test Kit
- Gram - Negative Bacteria Identification Test Kit
- Lactobacillus Identification Kit
- Listeria Identification Kit
- Neisseria Identification Kit
- Rapid Gram Screening Kit
- Salmonella Identification Kit
- Staph Identification Kit
- Strep Identification Kit
- Vibrio Identification Kit

**Presumptive Identification Tests**


- Rapid PYR Test Kit
- Rapid Presumptive Identification Test Kit for E.coli
- Rapid Test Kit for E.coli
- Test Kit for Esculin Hydrolysis Test Kit for MR-VP
- Rapid Test Kit for Phenylalanine Deamination
- Rapid Test Kit for Tryptophan Deamination

**Water Quality Testing Kits**

- PA Coliform Test Kit
- Coliform Test Kit
- H<sub>2</sub>S Test Kit
- H<sub>2</sub>S Medium Kit

**Chromogenic Media**

- Michrom™ Candida Differential Agar
- Michrom™ Coliform Agar
- Michrom™ Hydrocheck
- Michrom™ E. coli Agar
- Michrom™ Enterococcus Agar
- Michrom™ UTI Agar
- Michrom™ Vibrio Agar
- Michrom™ VRE Agar Base



**Blood Culturing System**

- BHI - Supplemented with 0.05% SPS
- Colorcult® Aerobic / Anaerobic / Paediatric / Selective Yeast and Fungi Culture Vial
- Colorcult®-System Pack (Aerobic / Paediatric Culture Vial)
- Glucose Broth - Supplemented with 0.05% SPS
- Soyabean Casein Digest Broth - Supplemented with 0.05%



**Micropro® System**

- Micropro® - BMD
- Micropro® - BCS (Broth Culture System)
- Micropro® - AST
- Micropro® - MIC
- Micropro® - ID



**Identification Of Bacteria**

**Easybact®**

**Biogram®**

Antimicrobial Susceptibility Discs as per CLSI / Customized Antimicrobial Susceptibility Discs




**Ready To Use Pouch Media**

- Instaprep® Nutrient Agar
- Instaprep® MacConkey Agar
- Instaprep® Mueller Hinton Agar
- Instaprep® Sabouraud Dextrose Agar
- Instaprep® Salmonella Shigella Agar



**Culture Media for Transporting Clinical Specimens**

- Amies with Charcoal Transport Medium
- Stuart Transport Medium
- Transport Swab with Cary Blair Medium



**Culture Media Recommended for Fungi, Yeast and Moulds**

- Corn Meal Agar
- Czapek Dox Agar
- Potato Dextrose Agar
- Sabouraud Chloramphenicol Agar
- Sabouraud Cycloheximide Chloramphenicol Agar
- Sabouraud Dextrose Agar

An ISO 13485 Certified Company  
**Microexpress®**

Division of

**T TULIP DIAGNOSTICS (P) LTD.**

Plot No. S-124, S-125, S-126, Utility Plot No. VIII, Phase III-B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403 202, INDIA.

Email: mex.queries@tulipgroup.com; Website: www.microexpress.in