


Microxpress[®]

Confident Results in 18 hours !




Growth Promotion

**Recovery
On XLD Agar**



Rappaport Vassiliadis Salmonella Enrichment Broth (Harmonized)

Selective enrichment medium for *Salmonella* species from pharmaceutical products in accordance with microbial testing by harmonized methodology of USP/EP/JP/IP.




Features	Benefits
1. Product formulated to perform as per international standards	<ul style="list-style-type: none">• Turbidity observed in 18 hours after inoculation of standard test strains.• Recovery on XLD Agar within 18 hours.
2. Wide storage temperature of 15°C-25°C	<ul style="list-style-type: none">• Better in use stability of product at user end.
3. Uniquely designed glass vials	<ul style="list-style-type: none">• Clear transparent material ensures better visualization of results.
4. Compact pack size of 25 x 10 mL vials	<ul style="list-style-type: none">• Ideally suited for small & big customers.• Occupies less storage space.
5. Ready to use product	<ul style="list-style-type: none">• Minimizes usage of equipment, manpower.

CAT No. : 203180240010 Pack Size : 25 x 10 mL Shelf Life : 1 Year

Bibliography

1. United State Pharmacopeia, 2019, Microbiological Examination of Nonsterile Products: Test for Specified Microorganisms <62>.
2. Data on file: Microxpress® , A Division of Tulip Diagnostics (P) Ltd.



An ISO 13485
Certified Company