

Steristik®

Gamma-Irradiated sterile cotton bud with Polypropylene stick.

Application

Steristik® cotton swabs are used for the collection of samples, surface sampling and may also be used for preparation of smears.

Technical Data Sheet

Parameters	UOM	Standard
Length	cm	13.5 ± 0.2
Outer diameter of polypropylene stick	cm	0.26 ± 0.02
Weight of the swab	g	0.35 ± 0.05
Colour		White

MOC

Stick - Polypropylene
Swab - Cotton

Note – Gamma-Irradiation Certificate provided on request.

Reference

Data on file: Microexpress®, A Division of Tulip Diagnostics (P) Ltd.

Product Presentation:

Cat No.	Product	Pack Size
208191270100	Steristik®	100 Nos.

Disclaimer

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.
