

Weiss Buffer**Intended Use**

Buffer for staining.

Summary & Principle

Weiss buffer solutions are used in hematology, for the dilution of staining solutions (e.g., Giemsa stock solution) and for the rinsing steps in the staining process.

Appearance

Clear, colourless solution

Storage and Stability

Store at 15-30°C away from bright light. Stability of the Weiss Buffer is as per expiry date mentioned on label.

Warranty

This product is designed to perform as described on the label and package insert. The manufacturer disclaims any implied warranty of use and sale for any other purpose.

Reference

1.Data on file: Microxpress®, A Division of Tulip Diagnostics (P) Ltd.

Product Presentation**Cat.No.**

207230070010

Product Description

Weiss Buffer

Pack Size

10mL

Disclaimer

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.
